

TRYPTONE BILE AGAR

A differential medium for the enumeration of *Escherichia coli* with DPM method.

Dehydrated media	
Code number:	500 g: TBA20500, 5 kg: TBA25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	7,0 – 7,4

Direction: Suspend **37 g** in one litre of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes.

Prepared media	
Bottled media:	100 ml: TBA30100, 500 ml: TBA30500
Plated media:	55 mm: TBA50055, 90 mm: TBA50090
Colour:	Yellowish
pH (25 °C):	7,1 – 7,3

Direction: Dispense the melted bottled media aseptically into sterile Petri-dishes. Media in Petri-dishes are ready to use.

FORMULA in g/l

Tryptone	20,5
Bile salts	1,5
Agar	15,0

Note: The typical formula can be adjusted to obtain optimal performance.

Storage conditions: Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled media protected from light at room temperature. Store the plated media protected from light at 2-8 °C. Use before the expiry date on the label.

Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Escherichia coli</i> ATCC 25922		Good	
<i>Staphylococcus aureus</i> ATCC 29213		Inhibited	

References: Anderson and Baird-Parker (1975) J. Appl. Bact. 39: 111.

For professional use only!