

ZIPPER BAG WITH WRITE-ON LABEL, STERILE

Resealable PE bag for transporting and storing materials

Name:	Zipper bag with write-on label, sterile
	for laboratory use only
Code number:	FZB93003

Specification:	
Packaging:	100 pcs/bag 3000 pcs/corrugated carton box
Material:	LDPE (low density polyethylene) bag with slide fastener
Colour:	Water clear
Size:	120 x 170 mm
Sterility:	Irradiated
Nominal treatment dose:	15 kGy
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.