

INOCULATING LOOP

Sterile loop with needle for the inoculation of micro-organism.

Name:	Inoculating loop, 10 µl, sterile
	for laboratory use only
Code number:	SLO91002

Specification:	
Packaging:	20 pcs/blister 1.000 pcs/unite 10.000 pcs/case
Packaging material:	Paper/foil blister
Colour:	Yellow
Length:	200 mm
Loop volume:	10 µl
Material:	Polystyrene
Sterility:	Irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.