

SPREADER

Sterile T-shape spreader for spreading the micro-organism on work surface.

Name:	Spreader
	for laboratory use only
Code number:	SPR90005

Specification:	
Packaging:	5 pcs/pouche 200 pcs/box
Packaging material:	BOPP foil
Colour:	Green
Material:	Polystyrene
Sterility:	Irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.