

SAMPLING SPOON

Plastic spoon for collecting samples.

Name:	Sampling spoon, individually wrapped, sterile
	for laboratory use only
Code number:	SSP90001

Specification:	
Packaging:	1.200 pcs/case 20 pcs/unit
Colour:	White
Volume:	5 ml
Material:	Polypropylene
Sterility:	Gamma irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.