

## PETRI-DISH

Plastic, disposable Petri-dish for culture media used for cultivation of micro-organisms in microbiology labs.

<b>Name:</b>	<b>Petri-dish 35x12,5 mm, thin wall, aseptic, triple vents</b>
	for laboratory use only
Code number:	PET90135

<b>Specification:</b>	
Packaging:	10 pcs/bag 160 pcs/case
Packaging material:	PE bag
Colour:	Water clear
Size:	35x12,5 mm
Wall thickness:	Thin
Vent:	Triple vented
Material:	Toxin free water clear polystyrene
Sterility:	Produced under aseptic conditions
SAL:	10 <sup>-2</sup>
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.