

PETRI-DISH

Plastic, disposable Petri-dish for culture media used for cultivation of micro-organisms in microbiology labs.

Name:	Petri-dish 55x14,6 mm, thin wall, aseptic, triple vents
	for laboratory use only
Code number:	PET90155

Specification:	
Packaging:	20 pcs/bag 1.420 pcs/case
Packaging material:	PE bag
Colour:	Water clear
Size:	55x14,6 mm
Wall thickness:	Thin
Vent:	Triple vented
Material:	Toxin free water clear polystyrene
Sterility:	Produced under aseptic conditions
SAL:	10 ⁻²
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.