

## PETRI-DISH

Plastic, disposable Petri-dish for culture media used for cultivation of micro-organisms in microbiology labs.

<b>Name:</b>	<b>Petri-dish 65 x 14 mm with grip lid, aseptic</b>
	for laboratory use only
Code number:	PET90165

<b>Specification:</b>	
Packaging:	20 pcs/bag, 720 pcs/box
Packaging material:	PE bag
Material - product:	Toxin free crystal-clear polystyrene
Material - collecting bag:	Low density polyethylene
Material - collecting box:	Corrugated carton
Colour:	Water clear
Size:	65x14 mm
Wall thickness:	Thin
Grip:	Yes
Material:	Toxin free water clear polystyrene
Heat resistance:	Not autoclavable
Sterility:	Manufactured aseptically
SAL:	10 <sup>-2</sup>

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.