

PETRI-DISH

Plastic, disposable Petri-dish for culture media used for cultivation of micro-organisms in microbiology labs.

Name:	Petri-dish 90x14,2 mm, thick wall, triple vents, 3 sectors, aseptic
	for laboratory use only
Code number:	PET90190-3S

Specification:	
Packaging:	25 pcs/bag 600 pcs/case
Packaging material:	PE bag
Colour:	Water clear
Size:	90x14,2 mm
Wall thickness:	Thick
Vent:	Triple vented
Sectors:	Three
Material:	Toxin free water clear polystyrene
Sterility:	Produced under aseptic conditions
SAL:	10 ⁻²
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.