

PETRI-DISH

Plastic, disposable Petri-dish for culture media used for cultivation of micro-organisms in microbiology labs.

Name:	Contact (Rodac) Petri-dish 65x14 mm, without grip lid, irradiated
	for laboratory use only
Code number:	PET90465

Specification:	
Packaging:	20 pcs/PE bag 720 pcs/corrugated carton box
Material:	Toxin free water clear polystyrene
Colour:	Water clear
Size:	65x14 mm
Wall thickness:	Thin
Grip:	No
Sterility:	Irradiated
Nominal treatment dose:	10 kGy
SAL:	10 ⁻⁵
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.