

PIPETTE TIP

Micropipette tips, made from high quality polypropylene. Tips are manufactured with the highest precision to ensure accurate and repeatable results in dispensing. The aerosol-resistant filter eliminates the air-to-liquid interface between the sample and the pipette.

Name:	Filter pipette tip 200 µl (Gilson type)
	for laboratory use only
Code number:	PIT91808

Specification:	
Packaging:	1.000 pcs/bag 50.000 pcs/case
Packaging material:	PE bag
Colour:	Yellow
Material, tip:	Polypropylene
Material, filter:	3-dimensional cross-linked sintered polyethylene
Heat resistance:	Not autoclavable
Suitable pipettes:	Gilson, Nichiryo, Finnpiquette, Biohit etc.
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.