

## PASTEUR PIPETTE, 1 ml, STERILE

The plastic Pasteur pipette is the more convenient alternative of the glass Pasteur pipettes. The product is disposable, chemically inert, non-toxic and significantly reduces the possibility of the contamination between the samples.

<b>Name:</b>	<b>Pasteur pipette, 1 ml, sterile, individually wrapped</b>
	for laboratory use only
<b>Code number:</b>	PPI92001

<b>Specification:</b>	
<b>Packaging:</b>	100 pcs/unit 2000 pcs/case
<b>Material:</b>	LDPE (low density polyethylene)
<b>Colour:</b>	Translucent
<b>Nominal fill volume:</b>	≈ 1 ml
<b>Volume of drops:</b>	≈ 50 µl
<b>Sterility:</b>	Irradiated
<b>SAL:</b>	10 <sup>-3</sup>
<b>Note:</b>	-

**Storage conditions:** Store in a dry place protected from light at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.