

## SEROLOGICAL PIPETTE

Disposable polystyrene serological pipettes aspirate and release millilitre liquid amounts using the blow-out method.

<b>Serological pipette, 1 ml, individually wrapped, sterile</b>	
	for laboratory use only
Code number:	SP-1-C

<b>Specification:</b>	
Packaging:	Individually wrapped - 2000 pcs/bag
Material - bag:	BOPP foil
Material - box:	Corrugated carton
Material:	Glass clear polystyrene
Length:	271 mm
Colour code:	Yellow
Nominal volume:	1,00 ml
Graduation:	0,01 ml
Heat resistance:	Not autoclavable
Sterility:	Irradiated
Nominal treatment dose:	20 kGy
SAL:	10 <sup>-5</sup>

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.