

SEROLOGICAL PIPETTE

Disposable polystyrene serological pipettes aspirate and release millilitre liquid amounts using the blow-out method.

Serological pipette, 2 ml, individually wrapped, sterile	
	for laboratory use only
Code number:	SP-2-C

Specification:	
Packaging:	Individually wrapped - 2000 pcs/bag
Material - bag:	BOPP foil
Material - box:	Corrugated carton
Material:	Glass clear polystyrene
Length:	271 mm
Colour code:	Green
Nominal volume:	2,00 ml
Graduation:	0,02 ml
Heat resistance:	Not autoclavable
Sterility:	Irradiated
Nominal treatment dose:	20 kGy
SAL:	10 ⁻⁵

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.