

COTTON TIPPED APPLICATOR

The cotton tipped applicator is a dry swab for many applications in microbiology labs, in environmental work, in veterinary labs, etc.

Name:	Cotton tipped applicator, wooden, sterile, individually wrapped
	for laboratory use only
Code number:	CTA90003

Specification:	
Packaging:	Individually wrapped 3000 pcs/case
Packaging material:	Paper
Material of stick:	Wooden
Dimensions of stick:	150 x 2,5 mm
Material of tip:	60% viscose + 40% pure cotton
Diameter of tip:	≈ 5 mm
Sterility:	Irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.