

## COTTON TIPPED APPLICATOR

The cotton tipped applicator is a dry swab for many applications in microbiology labs, in environmental work, in veterinary labs, etc.

<b>Name:</b>	<b>Cotton tipped applicator, wooden</b>
	for laboratory use only
Code number:	CTA90005

<b>Specification:</b>	
Packaging:	100 pcs/unit 3000 pcs/case
Packaging material:	PE bag
Material of stick:	Wooden
Dimensions of stick:	150 x 2 mm
Material of tip:	60% viscose + 40% pure cotton
Diameter of tip:	≈ 5 mm
Sterility:	Non-sterile
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.