

COTTON TIPPED APPLICATOR

The cotton tipped applicator is a dry swab for many applications in microbiology labs, in environmental work, in veterinary labs, etc.

Name:	Cotton tipped applicator, wire, sterile
	for laboratory use only
Code number:	CTA90009

Specification:	
Packaging:	Individually wrapped 3000 pcs/case
Packaging material:	Paper
Material of stick:	Aluminium
Dimensions of stick:	150 x 1 mm
Material of tip:	Cotton
Diameter of tip:	≈ 2 mm
Sterility:	Irradiated
Nominal treatment dose:	20 kGy
SAL:	10 ⁻⁵
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.