

COTTON TIPPED APPLICATOR

The cotton tipped applicator is a dry swab for many applications in microbiology labs, in environmental work, in veterinary labs, etc.

Name:	Cotton tipped applicator, wire, non sterile
	for laboratory use only
Code number:	CTA90010

Specification:	
Packaging:	100 pcs/unit 3000 pcs/case
Packaging material:	PE bag
Material of stick:	Aluminium
Dimensions of stick:	150 x 1 mm
Material of tip:	60% viscose + 40% pure cotton
Diameter of tip:	≈ 3 mm
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.