

SWAB APPLICATOR

Sterile dry swab in tube with cap, for taking different samples.

Name:	Swab applicator in tube, wooden, sterile
	for laboratory use only
Code number:	SWA90001

Specification:	
Packaging:	20 pcs/bag 2500 pcs/case
Packaging material:	PE bag with slide fastener
Material of tube and cap:	Polypropylene homopolymer
Material of stick:	Wooden
Dimensions of stick:	150 x 2,5 mm
Material of tip:	60% viscose + 40% pure cotton
Diameter of tip:	≈ 5 mm
Sterility:	Irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.