

SWAB APPLICATOR, IN TUBE, PLASTIC

Sterile dry swab in tube with cap, for taking different samples.

Name:	Swab applicator in tube, plastic, sterile
	for laboratory use only
Code number:	SWA90006

Specification:	
Packaging:	20 pcs/PE bag with slide fastener 2.500 pcs/corrugated carton box
Material of tube and cap:	Polypropylene
Material of stick:	High impact polystyrene
Material of tip:	60% viscose + 40% pure cotton
Dimensions of stick:	150 x 2,5 mm
Length of tip:	15 mm
Diameter of tip:	≈ 5 mm
Sterility:	Irradiated
Nominal treatment dose:	20 kGy
SAL:	10 ⁻⁵
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.