

## FAECAL CONTAINER

Plastic, disposable container, 3 parts (container + spoon + screw cap) for sampling, transport and storage of faeces samples.

<b>Name:</b>	<b>Faecal container, 3 parts</b>
	for laboratory use only
Code number:	FEC90002

<b>Specification:</b>	
Packaging:	1300 pcs/case
Colour of container:	Opalescent
Colour of cap:	White
Volume:	≈ 8 ml
Material:	Polypropylene – homopolymer
Heat resistance:	Autoclavable
Sterility:	Non sterile
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.