

CULTURE MEDIA TUBE WITH SPECIAL CUP

Disposable glass medium tube for sterilizing and storing culture media.

Name:	Culture Media Tube with Special Cup, 12 x 100 mm, in box
	for laboratory use only
Code number:	GTU91001

Specification:	
Packaging:	50 pcs/micro corrugated carton box 1800 pcs/cardboard box
Material of tube:	Glass
Colour of tube:	Water clear
Sizes:	100 x 12 mm
Volume:	≈ 10 ml
Heat resistance (tube and cup):	Autoclavable
Sterility:	Non-sterile

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.