

CULTURE MEDIA TUBE, 15 x 100 mm

Disposable glass medium tube for sterilizing and storing culture media.

Name:	Culture Media Tube, 15 x 100 mm
	for laboratory use only
Code number:	GTU91005

Specification:	
Packaging:	Bulk - 1500 pcs/cardboard box
Material of tube:	Glass
Colour of tube:	Water clear
Sizes:	100 x 15 mm
Volume:	≈ 13 ml
Sterility:	Non-sterile

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.