

CULTURE MEDIUM TUBE WITH SPECIAL CAP

Plastic, disposable tube, with special cap for sterilizing and storing culture media.

Name:	Culture medium tube with special cap 100x12 mm
	for laboratory use only
Code number:	TUB91001

Specification	
Packaging:	2000 pcs/corrugated box
Material:	High quality polypropylene random copolymer
Colour of tube:	Slightly opalescent
Colour of cap:	Dark blue
Size:	100 x 12 mm
Volume:	≈ 7 ml
Heat resistance:	Autoclavable
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.