

CULTURE MEDIUM TUBE WITH SPECIAL CAP

Plastic, disposable, sterile tube, with special cap for sterilizing and storing culture media.

Name:	Culture medium tube with special cap 100 x 12 mm, sterile
	for laboratory use only
Code number:	TUB92001

Specification	
Packaging:	100 pcs/bag, 3000 pcs/corrugated box
Material:	High quality polypropylene random copolymer
Colour of tube:	Slightly opalescent
Colour of cap:	Red
Size:	100 x 12 mm
Volume:	≈ 7 ml
Heat resistance:	Autoclavable
Sterility:	Irradiated
Nominal treatment dose:	10 kGy
SAL:	10 ⁻⁵
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.