

CULTURE MEDIUM TUBE WITH SPECIAL CAP

Plastic, disposable, sterile tube, with special cap for sterilizing and storing culture media.

| | |
|---------------------|--|
| Name: | Culture medium tube with special cap 150 x 16 mm, sterile |
| | for laboratory use only |
| Code number: | TUB92003 |

| Specification | |
|-------------------------|---|
| Packaging: | 25 pcs/bag, 1200 pcs/corrugated box |
| Material: | High quality polypropylene random copolymer |
| Colour of tube: | Slightly opalescent |
| Colour of cap: | Red |
| Size: | 150 x 16 mm |
| Volume: | ≈ 22 ml |
| Heat resistance: | Autoclavable |
| Sterility: | Irradiated |
| Nominal treatment dose: | 10 kGy |
| SAL: | 10 ⁻⁵ |
| Note: | - |

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.