

## EPPENDORF TUBE

Plastic, disposable, Eppendorf type, conical tube for centrifuge, transport and storage of biological samples.

<b>Name:</b>	<b>Eppendorf tube, 1,5 ml</b>
	for laboratory use only
Code number:	TUB92208

<b>Specification:</b>	
Packaging:	1.000 pcs/ PE bag; 10.000 pcs/corrugated box
Colour of container:	Slightly opalescent
Volume:	1,5 ml
Material:	High quality polypropylene - random copolymer
Heat resistance:	Autoclavable
Sterility:	Non sterile
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.