

CRYOTUBE

The vials, which are manufactured from high quality polypropylene, are suitable for storage of biological specimens at extremely low temperature.

Name:	Cryotube with write-on (frosted) area and with hub and gasket, 1,8 ml, sterile
	for laboratory use only
Code number:	TUB93001

Specification:	
Packaging:	500 pcs/bag 10.000 pcs/case
Packaging material:	PE bag
Colour:	Slightly opalescent
Volume:	1,8 ml
Material:	Special grade polypropylene
Heat resistance:	Autoclavable
Sterility:	Irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.