

## CRYOTUBE

The vials, which are manufactured from high quality polypropylene, are suitable for storage of biological specimens at extremely low temperature.

<b>Name:</b>	<b>Cryotube with write-on (frosted) area and with hub and gasket, 5 ml, sterile</b>
	for laboratory use only
Code number:	TUB93002

<b>Specification:</b>	
Packaging:	200 pcs/bag 3200 pcs/case
Packaging material:	PE bag
Colour:	Slightly opalescent
Volume:	5 ml
Material:	Special grade polypropylene
Heat resistance:	Autoclavable
Sterility:	Irradiated
SAL:	10 <sup>-3</sup>
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.