

URINE CONTAINER, 120 ML, STERILE

Transparent universal container 120 ml with inserted screw cap, sterile

Name:	Urine Container 120 ml, sterile
	for laboratory use only
Code:	URC90006

Specification:	
Packaging	300 pcs/box - inserted
Colour of cap:	Light blue
Colour of container:	Opalescent
Volume:	120 ml
Material	High quality, polypropylene random kopolimer
Sterility	Irradiated
Nominal treatment dose:	10 kGy
SAL	10 ⁻⁵

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.