

**URINE CONTAINER 60 ML,
SCREW CAP, STERILE**

Transparent universal container 60 ml with separated screw cap

Name:	Urine Container 60 ml, screw cap, sterile
	for laboratory use only
Code:	URC90008

Specification:	
Packiging	600 pcs/box
Colour of container:	Transparent
Volume:	60 ml
Material	High quality polypropylene
Sterility	Irradiated
SAL	10 ⁻³

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.