

URINE CONTAINER 120 ML, INDIVIDUALLY WRAPPED, STERILE

Transparent universal container 120 ml with inserted screw cap, individually wrapped, sterile

Name:	Urine Container 120 ml, individually wrapped, sterile
	for laboratory use only
Code:	URC90009

Specification:	
Packaging	Individually wrapped 250 pcs/box
Colour of cap:	Red
Colour of container:	Transparent
Volume:	120 ml
Material	High quality, water clear polypropylene
Sterility	Irradiated
Nominal treatment dose:	20 kGy
SAL	10 ⁻⁵

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.