

BRILLIANT GREEN

Brilliant green indicator for microbiological purposes.

Code number:	5 g: BRG10005, 25 g: BRG10025	
CAS: 633-03-4	Formula: $C_{27}H_{34}N_2O_4S$	Molecular Weight: 482,65

PHYSICO-CHEMICAL CHARACTERISTICS

Parameter	Specification
Appearance:	Dark green crystalline powder
Dye content [%]:	> 95
Absorption maximum in 50% ethanol [nm]:	628 – 632

Storage conditions: Store the product protected from light.