

BROMOCRESOL PURPLE

Bromocresol purple indicator for microbiological purposes.

Code number:	25 g: BCP10025, 100 g: BCP10100	
CAS: 143-74-8	Formula: C ₁₉ H ₁₄ O ₅ S	Molecular Weight: 540,24

PHYSICO-CHEMICAL CHARACTERISTICS

Parameter	Specification
Appearance:	Pink coloured powder
Visual transition range [from yellow to red]:	pH 5,2 - 6,8
Loss on drying [%]:	< 1,0

Storage conditions: Store the product protected from light.