

## KANAMYCIN SULPHATE

Kanamycin for microbiological purposes.

Code number:	500 g: KAN10500, 1000 g: KAN11000
Colour:	White
Appearance:	Fine powder, easily soluble in water

### Physico-chemical characteristics

Parameter	Specification
Assay of kanamycin	> 670 IU/mg
pH	5,5 - 7,5
Specific optical rotation	+103° to +115°
Loss on drying	< 5 %
Sulphated ash	< 0,5 %
Heavy metal	< 0,001 %

**Storage conditions:** Store the product in a dry place at room temperature.

**Antibiotic for in vitro diagnostic use – for professional use only!**