

## SAFRANINE

Safranine for microbiological purposes.

Code number:	5 g: SFR10025, 25 g: SFR10100	
CAS: 477-73-6	Formula: C <sub>20</sub> H <sub>19</sub> N <sub>4</sub> Cl	Molecular Weight: 350,84

### PHYSICO-CHEMICAL CHARACTERISTICS

Parameter	Specification
Appearance:	Dark greenish powder
Specific absorptivity [%/cm]:	> 875
Absorption maximum in 50% ethanol [nm]:	530 - 534
Loss on drying [%]:	< 15

**Storage conditions:** Store the product protected from light.