

## LISTERIA SELECTIVE SUPPLEMENT, UVM I

FOR 500 ml OF LISTERIA ENRICHMENT BROTH, UVM I

Freeze-dried mixture for the isolation of *Listeria monocytogenes*.

Description	
Code number:	4 ml: LU180004
Colour:	Yellowish
Appearance:	Homogeneous lyophilisate

**Direction:** Suspend **27,5 g of Listeria Enrichment Broth Base, UVM-Fraser (LEF20500)** in 500 ml of distilled water and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial of Listeria Selective Supplement, UVM I (LU180004)** reconstituted with 4 ml of 1:1 mixture of ethanol and sterile distilled water. Mix well. Dispense aseptically into sterile final containers.

### FORMULA [mg/vial]

Nalidixic acid	10
Acriflavine	6

**Storage conditions:** Protected from light, at 2-8 °C.

#### Warning!

Use before the expiry date on the label.

**For professional use only!**