

## TWEEN 80 SUPPLEMENT

TWEEN 80 (Polyoxyethylene sorbitan monooleate; Polysorbate 80) for the preparation of some media.

Description	
Code number:	100 ml: TWS80100, 500 ml: TWS80500 2500 ml: TWS82500
Colour:	Colourless or slightly yellowish
Appearance:	Thick, transparent liquid. A little precipitation may occur but not detrimental.

**Direction:** Various. See the product information of the relevant media bases.

### FORMULA

TWS80100	TWEEN 80 (bacteriological grade) - 100 ml
TWS80500	TWEEN 80 (bacteriological grade) - 500 ml
TWS82500	TWEEN 80 (bacteriological grade) - 2500 ml

**Storage conditions:** Protected from light, at room temperature.

**Warning!**

Use before the expiry date on the label.

**For professional use only!**